



Certificate of Laboratory Testing

Issue No.	R20191220-0007	Receipt No.	190200538-001
Test Completion Date	2019-12-20	Date of Receipt	2019-11-05
Product Name	:dc Microparticle Stopper mask KF94(Gray mask with Yellow ear loops, L size)		
Manufacturing report number		Batch Number	Dccx02
Product Type, Material, Classification	Mask		
Date of Manufacturing	2019-09-12		
Expiration Date			
Client	Company Name	Dacian Technology Material Co.,Ltd	
	Name	Fangling hung	
	Address	No8, Lan 2438, JingQuo RD, Dajai, 437, Taichung city, Taiwan (Tel. No.:886-4-26810387 Fax No.: E-mail:fanny@dctpro.com)	
Manufacturer	Name	Dacian Technology Material Co.,Ltd	
	Country	Taiwan	
	Address	No8, Lan 2438, JingQuo RD, Dajai, 437, Taichung city, Taiwan	
Test Purpose	Import inspection		
Remarks			

Test Items and Results

Test Items	Specifications	Test Results	Judgement	Remarks
Visual	PASS/FAIL	PASS	Suitability	
Shape	PASS/FAIL	PASS	Suitability	
Leakage rate test	11.0% and below	PASS	Suitability	
Breathing resistance test	70 Pa and below	PASS	Suitability	
Particle filter penetration test (Sodium chloride test method)	94 % and above	PASS	Suitability	
Particle filter penetration test (Paraffin oil test method)	94 % and above	PASS	Suitability	



Test Items	Specifications	Test Results	Judgement	Remarks
Purity(Acid and alkali)	PASS/FAIL	PASS	Suitability	
Purity(Identification of Fluorescent Brightening Agents)	PASS/FAIL	PASS	Suitability	
Purity(Pigment)	PASS/FAIL	PASS	Suitability	
Purity(Formaldehyde)	PASS/FAIL	PASS	Suitability	
Tensil strength of hairband joint	10N and above	PASS	Suitability	

Final Judgment : Suitability

Testing Personnel: Kwon Dae Hyuk, Lee So Yeon

Testing manager : Kong Yong Won

- ※ The judgment above shall be only subject to the requested test items.
- ※ In the case of insufficient space on the page, test items and test results may be written in the appendix.
- ※ When advertising test results or indicating test results on containers or packages, the entire content of the test report shall be indicated

The above testing and inspection report shall be issued in accordance with Article 11(2) of the Act on Testing and inspection of Food and Drugs and Article 12(4)1 of the Enforcement Regulation of the Act

20 / 12 / 2019

